TIB D-1003

TECHNICAL INFORMATION BULLETIN

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Connection of VarioLine devices to medifon 2010 / 2000 communications terminal

This Bulletin sets out the possible options for using Medical 800 VarioLine devices in combination with the KT2010 / KT2000 communications terminal.

Please refer to the table below for the various combination options and their features.

VL devices with base module	Possible patient's handsets	KT2010/2000	KT2010 function	KT2000 function
125.010x Base module with button	none	yes / yes	Call triggered according to voltage coding of button	Triggers a 24 V call, button must be coded to 24 V
125.050x Base module with socket	125.5600 125.5610 125.5620	yes / yes	Call triggered according to voltage coding of button, speech possible via KT	Triggers a 24 V call, button must be coded to 24 V
125.051x Base module with socket, 15p, VDE	125.5600 125.5610 125.5620	yes / yes	Call triggered according to voltage coding of button	Triggers a 24 V call, button must be coded to 24 V
125.055x Base module with socket, 15p, ribbon cable with 125.080x VL Adapter with radio, 15p	125.5600 125.5610 125.5620 125.5720	yes / yes	Call triggered according to voltage coding of button, speech possible via KT, speech via patient's handset only possible with additional PCB*)	Triggers a 24 V call, button must be coded to 24 V, speech possible via KT2000
125.060x Pull switch, 125.061x Pneumatic switch	none	yes / yes	Triggers a 24 V call, button must be coded to 24 V	Triggers a 24 V call, button must be coded to 24 V
125.0701 Base module with socket for MSU, 15p, VDE with 125.080x VL Adapter with radio, 15p	125.5600 125.5610 125.5620 125.5720	yes / yes	Call triggered according to voltage coding of button, speech possible via KT, speech via patient's handset only possible with additional PCB*)	Triggers a 24 V call, button must be coded to 24 V, speech possible via KT2000

*) Additional PCB

The VL Adapter with radio 125.080x does not currently support bedwise speech if connected to a KT 2010/KT2000. If this function is to be realised anyway, however, an additional PCB is required. This must be connected to the VL Adapter. The wiring diagram can be obtained from TSP Munich Care Communications Support on request.

Principle

All combination options not listed in the above table are not approved by TSP Munich.

Outlook

In the next redesign of the VL Adapter with radio 15p the additional circuitry currently on the additional PCB will be incorporated into the VL Adapter. We will inform you by bulletin as soon as this modified VL Adapter is available.

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Michael Koschek
09.02.2010

For further information, please contact Michael Koschek, Tyco Safety Products – Technical Support Engineer, Munich, Germany (mkoschek@tycoint.com).

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